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
Inhibitory RNaz z rodziny RNasin®

Table 1. Promega RNA Analysis Products That Are Compatible with RNasin® Ribonuclease Inhibitors.

Compatible Application	Product	Cat.#
RT-PCR	Access RT-PCR System ^(a,d)	A1250
	AccessQuick™ RT-PCR System ^(d)	A1702
First-Strand cDNA Synthesis	✓ Reverse Transcription System ^(a,b)	A3500
	✓ ImProm-II™ Reverse Transcription System ^(a,b)	A3800
First- and Second-Strand cDNA Synthesis	Universal RiboClone® cDNA Synthesis System ^(a,b,e)	C4360
	✓ Synthesis System ^(a,b,e)	C4360
Reverse Transcriptases	M-MLV RT	M1701
	M-MLV RT RNase H–	M5301
	M-MLV RT RNase H– Point Mutant	M3681
	AMV RT	M5101
	ImProm-II™ Reverse Transcriptase	A3802
RNA Polymerases	SP6 RNA Polymerase	P1081
	T3 RNA Polymerase	P2083
	T7 RNA Polymerase	P2075
In vitro Transcription Systems	SP6 Riboprobe® System ^(a,b)	P1420
	T3 Riboprobe® System ^(a,b)	P1430
	T7 Riboprobe® System ^(a,b)	P1440
Large-Scale RNA Production	✓ T7 RiboMAX™ Express RNAi System ^(a,b,f)	P1700
	✓ T7 RiboMAX™ Express Large Scale RNA Production System ^(a,b)	P1320
	✓ RiboMAX™ Large Scale RNA Production System—SP6 ^(a,b,g)	P1280
	✓ RiboMAX™ Large Scale RNA Production System—T7 ^(a,b,g,h)	P1300
In vitro Translation	Rabbit Reticulocyte Lysate ^(a,g,i)	L4960
	Wheat Germ Extract	L4380
Coupled Transcription/Translation	TnT® SP6 Coupled Reticulocyte Lysate System ^(a,g,i,j)	L4600
	TnT® T3 Coupled Reticulocyte Lysate System ^(a,g,i,j)	L4950
	TnT® T7 Coupled Reticulocyte Lysate System ^(a,g,i,j)	L4610
	TnT® SP6 Coupled Wheat Germ Extract System ^(a,g,i,j)	L4130
	TnT® T3 Coupled Wheat Germ Extract System ^(a,g,i,j)	L4120
Coupled Transcription/Translation Master Mix	✓ TnT® SP6 Quick Coupled Transcription/Translation System ^(a,b,g,i,j)	L2080
	✓ TnT® T7 Quick Coupled Transcription/Translation System ^(a,b,g,i,j)	L1170
	✓ TnT® T7 Quick for PCR DNA ^(a,b,i)	L5540

✓ Denotes that system includes Recombinant RNasin® Ribonuclease Inhibitor either as a separate tube or within the Master Mix.

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- **RNasin® Ribonuclease Inhibitors Natural (N2111, N2115)**

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- **RNasin® Plus RNase Inhibitor New Protein for High Temperature RNase Inhibition (N2611, N2615)**

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Table 1. Effect of Natural and Recombinant RNasin® Ribonuclease Inhibitors on Selected Nucleases and Polymerases.

Inhibits	Does Not Inhibit
RNase A	RNase T1
RNase B	S1 Nuclease
RNase C	RNase from <i>Aspergillus</i>
human placental RNase	RNase H
angiogenin	RNase ONE™ Ribonuclease
	Taq DNA Polymerase
	AMV or M-MLV Reverse Transcriptase
	SP6, T7, T3 RNA Polymerase

Product	Size	Cat.#
Recombinant RNasin® Ribonuclease Inhibitor ^(a,b)	2,500u	N2511
	10,000u	N2515
RNasin® Plus RNase Inhibitor ^(c)	2,500u	N2611
	10,000u	N2615

For Laboratory Use.